

■ SYRINGE PUMP MR-QD



Continuous flow from 750 µl/h up to 15 l/h
(depending on syringe volume)

- Max. pressure: 145 psi
- Fluids in contact with PTFE and Perlast®
- Syringes from 0,25 ml up to 25 ml
- Control unit, pressure and temperature sensor included
- Control unit can manage up to 4 pumps
- Self-priming
- Self-optimization pulsation mode
- Safety stop on pressure and temperature
- Various operation modes
- Suitable for Flow Chemistry Starter Set

